

## Rat PKCa (Protein Kinase C Alpha) Sandwich ELISA Kit Ready-To-Use Kit (STJE0017883)

STJE0017883

### GENERAL INFORMATION


<b>Product</b>	Kits
<b>Type</b>	
<b>Short</b>	This PKCa Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat
<b>Description</b>	tissue homogenates, cell lysates or other biological fluids..
<b>Applications</b>	ELISA
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Detection</b>	0.625-40ng/mL
<b>Sensitivity</b>	0.237ng/mL
<b>Detection</b>	0.625-40ng/mL
<b>Limit</b>	
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Instruction</b>	
<b>Note</b>	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

### TARGET INFORMATION

<b>Sample</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Types</b>	
<b>Gene</b>	<a href="#">Prkca</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">KPCA_RAT</a>



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081