

## Human SCT (Secretin) Competitive ELISA Kit Ready-To-Use Kit (STJE0018432)

STJE0018432

### GENERAL INFORMATION


|                          |                                                                                                                                                                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>           | Kits                                                                                                                                                                                                                                        |
| <b>Type</b>              |                                                                                                                                                                                                                                             |
| <b>Short Description</b> | This SCT Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.. |
| <b>Applications</b>      | ELISA                                                                                                                                                                                                                                       |
| <b>Reactivity</b>        | Human                                                                                                                                                                                                                                       |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | 23.4-1500pg/mL                                                                                                                                                                                                                                 |
| <b>Sensitivity</b>     | 7.25pg/mL                                                                                                                                                                                                                                      |
| <b>Detection Limit</b> | 23.4-1500pg/mL                                                                                                                                                                                                                                 |
| <b>Assay Time</b>      | 2.5 hrs                                                                                                                                                                                                                                        |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Instruction</b>     |                                                                                                                                                                                                                                                |
| <b>Note</b>            | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                                                                                        |

### TARGET INFORMATION

|                   |                                                                                                      |
|-------------------|------------------------------------------------------------------------------------------------------|
| <b>Sample</b>     | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Types</b>      |                                                                                                      |
| <b>Gene</b>       | SCT                                                                                                  |
| <b>Symbol</b>     |                                                                                                      |
| <b>Uniprot ID</b> | SECR_HUMAN                                                                                           |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081