

## Mouse TUBa4A (Tubulin Alpha 4A) Sandwich ELISA Kit Ready-To-Use Kit (STJE0019250)

STJE0019250

### GENERAL INFORMATION


|                          |                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                                                |
| <b>Short Description</b> | This TUBa4A Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse tissue homogenates, cell lysates or other biological fluids.. |
| <b>Applications</b>      | ELISA                                                                                                                                                                                               |
| <b>Reactivity</b>        | Mouse                                                                                                                                                                                               |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | 0.312-20ng/mL                                                                                                                                                                                                                                  |
| <b>Sensitivity</b>     | 0.114ng/mL                                                                                                                                                                                                                                     |
| <b>Detection Limit</b> | 0.312-20ng/mL                                                                                                                                                                                                                                  |
| <b>Assay Time</b>      | 3 hrs                                                                                                                                                                                                                                          |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Instruction</b>     |                                                                                                                                                                                                                                                |
| <b>Note</b>            | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                                                                                        |

### TARGET INFORMATION

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| <b>Sample Types</b> | tissue homogenates, cell lysates or other biological fluids. |
| <b>Gene</b>         | <a href="#">Tuba4a</a>                                       |
| <b>Symbol</b>       |                                                              |
| <b>Uniprot ID</b>   | <a href="#">TBA4A_MOUSE</a>                                  |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081