

## Mouse IL-1 beta (IL-1 beta) Sandwich ELISA Kit Kit (STJE0019551)

STJE0019551

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	Colorimetric
<b>Sensitivity</b>	7.8 pg/mL
<b>Detection Limit</b>	156-1000 pg/ml
<b>Assay Time</b>	4 hours
<b>Storage</b>	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
<b>Instruction Note</b>	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma and cell culture supernates
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	