

## Mouse IgG (Total) (IgG (Total)) Sandwich ELISA Kit Kit (STJE0019591)

STJE0019591

### GENERAL INFORMATION

|                          |       |
|--------------------------|-------|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> |       |
| <b>Applications</b>      | ELISA |
| <b>Reactivity</b>        | Mouse |

### PRODUCT PROPERTIES

|                         |                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>        | Colorimetric                                                                                                            |
| <b>Sensitivity</b>      | 0.08 ng/mL                                                                                                              |
| <b>Detection Limit</b>  | 156-100 ng/ml                                                                                                           |
| <b>Assay Time</b>       | 4 hours                                                                                                                 |
| <b>Storage</b>          | Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.                       |
| <b>Instruction Note</b> | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |

### TARGET INFORMATION

|                     |                                           |
|---------------------|-------------------------------------------|
| <b>Sample Types</b> | serum, plasma and cell culture supernates |
| <b>Gene Symbol</b>  |                                           |
| <b>Uniprot ID</b>   |                                           |