

Mouse Monocyte Chemotactic Protein 1 (Monocyte Chemotactic Protein 1) Sandwich ELISA Kit (STJE0019593)

STJE0019593

GENERAL INFORMATION

Product Type	Kits
Short Description	
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	1.5 pg/mL
Detection Limit	312-200 pg/ml
Assay Time	4 hours
Storage	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and cell culture supernates
Gene Symbol	
Uniprot ID	