

Human Lactoferrin (LTF) (Lactoferrin (LTF)) Sandwich ELISA Kit Kit (STJE0019609)

STJE0019609

GENERAL INFORMATION

Product Type	Kits
Short Description	
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.39 ng/mL
Detection Limit	078-50 ng/mL
Assay Time	4 hours
Storage	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma, cell culture supernatants, cell lysates, tissue homogenates and other biological fluids
Gene Symbol	
Uniprot ID	