

Sars-Cov-2 Nucleocapsid Protein IgG Antibody (Nucleocapsid Protein IgG Antibody) Indirect ELISA Kit Kit (STJE0019637)

STJE0019637

GENERAL INFORMATION

Product Type	Kits
Short Description	
Applications	ELISA
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	
Detection Limit	
Assay Time	3.5 hours
Storage	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma
Gene Symbol	
Uniprot ID	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081