

Spike S1 Protein Antibody (Spike S1 Protein Antibody) Sandwich ELISA Kit Kit (STJE0019639)

STJE0019639

GENERAL INFORMATION


Product Type Kits
Short Description
Applications
Reactivity

PRODUCT PROPERTIES

Detection Colorimetric
Sensitivity 182.1 pg/mL
Detection Limit 078-50 ng/mL
Assay Time 4 hours
Storage Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction
Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types serum, plasma, cell culture supernates and other biological fluids
Gene Symbol
Uniprot ID



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081