

## Mouse INS (Insulin) Sandwich ELISA Kit (STJE0019689)

STJE0019689

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This INS Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant and tissue homogenate.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	0.19ng/mL
<b>Detection Limit</b>	0.31-20ng/mL
<b>Assay Time</b>	3.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Instruction</b>	
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Cell culture supernatant and tissue homogenate
<b>Gene Symbol</b>	<a href="#">Ins2</a> <a href="#">Ins1</a>
<b>Uniprot ID</b>	<a href="#">INS2_MOUSE</a> <a href="#">INS1_MOUSE</a>