

Rat INS (Insulin) Sandwich ELISA Kit (STJE0019690)

STJE0019690

GENERAL INFORMATION

Product Type	Kits
Short Description	This INS Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant and tissue homogenate.
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.19ng/mL
Detection Limit	0.31-20ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Cell culture supernatant and tissue homogenate
Gene Symbol	Ins2 Ins1
Uniprot ID	INS2_RAT INS1_RAT