

Mouse TNF-Alpha (Tumor Necrosis Factor Alpha) Sandwich ELISA Kit (STJE0019716)

STJE0019716

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNF-Alpha Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	4.69pg/mL
Detection Limit	7.81-500pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Tnf
Uniprot ID	TNFA_MOUSE