

Human DBP (Vitamin D-binding protein) Sandwich ELISA Kit (STJE0019757)

STJE0019757

GENERAL INFORMATION

Product Type	Kits
Short Description	This DBP Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.09ng/mL
Detection Limit	0.16~10ng/mL
Assay Time	3.5h
Storage Instruction	An unopened kit can be stored at 2-8 °C for 1 week. If the kit is not supposed to be used within 1 week, store the items separately according to the following conditions once the kit is received.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	GC
Symbol	
Uniprot ID	VTDB_HUMAN