

Mouse DPP4 (Dipeptidyl Peptidase 4) Sandwich ELISA Kit (STJE0019797)

STJE0019797

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This DPP4 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids. |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | Colormetric |
| Sensitivity | 0.94ng/mL |
| Detection Limit | 1.56~100ng/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma and other biological fluids |
| Gene Symbol | Dpp4 |
| Uniprot ID | DPP4_MOUSE |