

Rat IL-1R2 (Interleukin 1 Receptor Type 2) Sandwich ELISA Kit (STJE0019802)

STJE0019802

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL-1R2 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	18.75pg/mL
Detection Limit	31.25~2000pg/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and other biological fluids
Gene Symbol	IL1r2
Uniprot ID	IL1R2_RAT