

Human MIP-1 Alpha (Macrophage Inflammatory Protein 1 Alpha) Sandwich ELISA Kit (STJE0019880)

STJE0019880

GENERAL INFORMATION

Product Type	Kits
Short Description	This MIP-1 Alpha Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	1.62 pg/mL
Detection Limit	3.13-200 pg/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and other biological fluids
Gene Symbol	CCL3
Uniprot ID	CCL3_HUMAN