

## Human TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) Sandwich ELISA Kit (STJE0019893)

STJE0019893

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This TNFSF14 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	2.67 pg/mL
<b>Detection Limit</b>	7.81-500 pg/mL
<b>Assay Time</b>	3.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Gene Symbol</b>	<a href="#">TNFSF14</a>
<b>Uniprot ID</b>	<a href="#">TNF14_HUMAN</a>