

Rat NGAL (Neutrophil Gelatinase Associated Lipocalin) Sandwich ELISA Kit (STJE0019945)

STJE0019945

GENERAL INFORMATION

Product Type	Kits
Description	This NGAL Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.84 pg/mL
Detection Limit	1.56-100 pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and other biological fluids
Gene Symbol	Lcn2
Uniprot ID	NGAL_RAT