

Rat IgG2b (Immunoglobulin G2b) Sandwich ELISA Kit (STJE0019950)

STJE0019950

GENERAL INFORMATION

| | |
|---------------------|---|
| Product Type | Kits |
| Description | This IgG2b Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma and other biological fluids. |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | Colormetric |
| Sensitivity | 65.63 pg/mL |
| Detection Limit | 109.38-7000 pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma and other biological fluids |
| Gene Symbol | |
| Uniprot ID | |