

Monkey IGF-1 (IGF-1) Sandwich ELISA Kit (STJE0020058)

STJE0020058

GENERAL INFORMATION

Product Type	Kits
Short Description	This IGF-1 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in monkey serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Monkey

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	34.56 pg/mL
Detection Limit	93.75-6000 pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	