

Porcine IL-1 Alpha (IL-1 Alpha) Sandwich ELISA Kit (STJE0020073)

STJE0020073

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL-1 Alpha Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in porcine serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Porcine

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	9.38 pg/mL
Detection Limit	15.63-1000 pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	