

## Pig GDF15 (GDF15) Sandwich ELISA Kit (STJE0020076)

STJE0020076

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This GDF15 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in pig serum, plasma and other biological fluids. |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Pig   |

### PRODUCT PROPERTIES

|                        |  |
|------------------------|--|
| <b>Detection</b>       | Colormetric  |
| <b>Sensitivity</b>     | 14.78 pg/mL  |
| <b>Detection Limit</b> | 39.06-2500 pg/mL   |
| <b>Assay Time</b>      | 3.5h   |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |  |
| <b>Note</b>            | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.  |

### TARGET INFORMATION

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  |   |
| <b>Uniprot ID</b>   |   |