

2019-nCoV 3C-like Proteinase protein (His-Tag) (STJP000001)

STJP000001

GENERAL INFORMATION

Product Type	Proteins
Short Description	2019-nCoV 3C-like Proteinase-His-Tag protein was developed from and has a target region of His-Tag. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	
Formulation	
Purification	
Dilution Range	>90% as determined by SDS-PAGE
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081