

2019-nCoV Nucleocapsid protein protein (His-Tag) (STJP000006)

STJP000006

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | 2019-nCoV Nucleocapsid protein-His-Tag protein was developed from and has a target region of His-Tag. For use in research applications. |
| Host/Source | |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | |
| Purification | |
| Dilution Range | >90% as determined by SDS-PAGE |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|--------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Sequence |