

2019-nCoV Spike protein fragment 2 protein (His-Tag) (STJP000009)

STJP000009

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | 2019-nCoV Spike protein fragment 2-His-Tag protein was developed from and has a target region of His-Tag. For use in research applications. |
| Host/Source | |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | |
| Purification | |
| Dilution Range | >90% as determined by SDS-PAGE |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|--------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Sequence |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081