

2019-nCoV Spike S1 protein (C-His) (STJP000023)

STJP000023

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. or Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
Purification
Dilution Range >90% by SDS-PAGE;> 95% by HPLC
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081