

## ACE2 protein (C-His&Avi) (STJP000026)

STJP000026

### GENERAL INFORMATION





**Product Type** Proteins  
**Short Description**  
**Host/Source** HEK293 cells

### PRODUCT PROPERTIES

**Concentration**  
**Formulation** Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. or Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.  
**Purification**  
**Dilution Range** > 95% by SDS-PAGE.  
**Storage Instruction** Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

**Gene ID** 59272  
**Gene Symbol** ACE2  
**Uniprot ID** ACE2\_HUMAN  
**Immunogen Sequence**



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081