

MANF protein (STJP000139)

STJP000139

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | MANF-protein was developed from. For use in research applications. |
| Host/Source | |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | 1 mg/ml |
| Formulation | PBS pH 7.4, concentration |
| Purification | Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media |
| Dilution Range | |
| Storage Instruction | Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 7873 |
| Gene Symbol | MANF |
| Uniprot ID | MANF_HUMAN |
| Immunogen Sequence | |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081