

2B4 protein (STJP000140)

STJP000140

GENERAL INFORMATION

Product Type	Proteins
Short Description	2B4-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	PBS pH 6.0
Purification	Purified by Ni-affinity chromatography, ion-exchange chromatography and desalting chromatography from serum-free CHO growth media
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
Uniprot ID	CD244_HUMAN
Immunogen	
Sequence	