

2B4 protein (STJP000141)

STJP000141

GENERAL INFORMATION

Product Type	Proteins
Short Description	2B4-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	PBS pH 6.0
Purification	Purified by Ni-affinity chromatography, ion-exchange chromatography and desalting chromatography from serum-free CHO growth media
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	18106
Gene Symbol	Cd244
Uniprot ID	CD244_MOUSE
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081