

GNDF protein (STJP000148)

STJP000148

GENERAL INFORMATION

Product Type Proteins
Short Description GDNF-protein was developed from. For use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 0.96 mg/ml
Formulation PBS pH 7.4
Purification Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media. Protein is sterile-filtrated through 0.22 µm filter
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [2668](#)
Gene Symbol [GDNF](#)
Uniprot ID [GDNF_HUMAN](#)
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081