

## GFRA-1 protein (STJP000149)

STJP000149

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	GFRA-1-protein was developed from. For use in research applications.
<b>Host/Source</b>	

### PRODUCT PROPERTIES

<b>Concentration</b>	0.9 mg/ml
<b>Formulation</b>	PBS pH 7.4
<b>Purification</b>	Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>	2674
<b>Gene Symbol</b>	GFRA1
<b>Uniprot ID</b>	GFRA1_HUMAN
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081