

GFRA-2 protein (STJP000150)

STJP000150

GENERAL INFORMATION

Product Type Proteins
Short Description GFRA-2-protein was developed from. For use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 0.7 mg/ml
Formulation PBS pH 7.4
Purification
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [2675](#)
Gene Symbol [GFRA2](#)
Uniprot ID [GFRA2_HUMAN](#)
Immunogen Sequence