

GFRA-2 protein (STJP000153)

STJP000153

GENERAL INFORMATION

Product Type	Proteins
Short Description	GFRA-2-protein was developed from. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	1.3 mg/ml
Formulation	PBS pH 7.4
Purification	Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	2676
Gene Symbol	GFRA3
Uniprot ID	GFRA3_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081