

ACE2 protein (hFc) (STJP000155)
STJP000155**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ACE2-hFc protein was developed from human and has a target region of hFc. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	1. PBS pH 7.4, 2. 20 mM Na-citrate, 150 mM NaCl; pH 6.0
Purification	Purified by Protein A affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. NA
Dilution Range	> 95%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
Uniprot ID	ACE2_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081