

ACE2 protein (His-Tag) (STJP000156)

STJP000156

GENERAL INFORMATION

Product Type	Proteins
Short Description	ACE2-His-Tag protein was developed from human and has a target region of His-Tag. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. NA
Dilution Range	> 80%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
Uniprot ID	ACE2_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081