

RNase 7 protein (STJP000162)

STJP000162

GENERAL INFORMATION

Product Type	Proteins
Short Description	RNase 7-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1.9 mg/ml
Formulation	10 mM Na-acetate, pH 5.0; 150 mM NaCl
Purification	Purified by affinity chromatography and gel-filtration from serum-free CHO growth media.
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84659
Gene Symbol	RNASE7
Uniprot ID	RNASE7_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081