

## Recombinant SARS-CoV-2 Spike S1 protein (STJP000164)

STJP000164

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant SARS-CoV-2 Spike S1-protein was developed from sars-cov-2. For use in research applications.
<b>Host/Source</b>	SARS-CoV-2

### PRODUCT PROPERTIES

<b>Concentration</b>	1 mg/ml
<b>Formulation</b>	PBS pH 7.4
<b>Purification</b>	NA Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. NA
<b>Dilution Range</b>	>90%
<b>Storage Instruction</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen Sequence</b>

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081