

Recombinant SARS-CoV-2 Spike S2 ECD-His protein (STJP000165)

STJP000165

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant SARS-CoV-2 Spike S2 ECD-His-protein was developed from sars-cov-2. For use in research applications.
Host/Source	SARS-CoV-2

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. NA
Dilution Range	> 90%
Storage	Store at -70°C upon receipt. It is strongly recommended to avoid repeated freeze-thaw cycles. Aliquot the protein, freeze quickly in liquid nitrogen and thaw only once for your experiments.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081