

Recombinant SARS-CoV-2 Trimeric Spike protein (STJP000166)

STJP000166

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant SARS-CoV-2 Trimeric Spike-protein was developed from. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated
Dilution Range	>95%
Storage Instruction	Store at -70°C. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081