

SARS-CoV-2 Trimeric Spike Indian B.1.617.1 protein (Kappa) (STJP000168)

STJP000168

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 Trimeric Spike Indian B.1.617.1-Kappa protein was developed from cho cells and has a target region of Kappa. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated NA
Dilution Range	>95%
Storage	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence