

SARS-CoV-2 Trimeric Spike UK VOC 202012/01 protein (Alpha) (STJP000177)

STJP000177

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 Trimeric Spike UK VOC 202012/01 -Alpha protein was developed from cho cells and has a target region of Alpha. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated
Dilution Range	>95%
Storage	Stored at +4 °C
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081