

SARS-CoV-2 Spike S1 UK VOC 202012/01 protein (Alpha) (STJP000182)

STJP000182

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 Spike S1 UK VOC 202012/01-Alpha protein was developed from and has a target region of Alpha. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter.
Dilution Range	>95%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081