

Human CD16b (FcγR3b) protein (Recombinant) {Biotin} (STJP000189)

STJP000189

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD16b (FcγR3b)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised 0.2 µg/ml filtered PBS solution, pH7.4 containing 5% trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	2215
Gene Symbol	FCGR3B
Uniprot ID	FCG3B_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081