

Human M-CSF protein (Recombinant) {Biotin} (STJP000195)

STJP000195

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human M-CSF-protein was developed from b.mori. For use in research applications.
Host/Source	B.mori

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered PBS solution pH7.0 containing % HSA, % mannitol.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
Instruction	loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	1435
Gene Symbol	CSF1
Uniprot ID	CSF1_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081