

Mouse VEGF 164AA protein (Recombinant) (STJP000241)

STJP000241

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse VEGF 164AA-protein was developed from. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration
Formulation
Purification
Dilution Range
Storage Instruction

TARGET INFORMATION

Gene ID	22339
Gene Symbol	Vegfa
Uniprot ID	VEGFA_MOUSE
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081