

Hepatitis C Virus Nonstructural protein 4 protein (Recombinant) (STJP000250)

STJP000250

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Hepatitis C Virus Nonstructural protein 4-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 M Tris solution pH7.5.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at -20°C for one month, or at -80°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081