

**MERS Nucleocapsid protein protein (Recombinant) (STJP000257)**

STJP000257

**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-MERS Nucleocapsid protein-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

**PRODUCT PROPERTIES**

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC NA
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

**TARGET INFORMATION**

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen Sequence**

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081