

2019 nCoV/COVID19 Spike protein S protein (Recombinant) (STJP000265)

STJP000265

GENERAL INFORMATION


Product Type Proteins
Short Description Recombinant-2019 nCoV/COVID19 Spike protein S-protein was developed from hek293. For use in research applications.
Host/Source HEK293

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from 0.2 µm filtered 20 mM Tris, 300 mM NaCl, 10 % glycerol, pH 8.0.
Purification
Dilution Range >95%, as determined by SDS-PAGE and HPLC NA
Storage Instruction The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



1µg loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081