

Mouse IL-13 alpha 1 receptor (IL13RA1) protein (Recombinant) (STJP000279)

STJP000279

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-13 alpha 1 receptor (IL13RA1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	16164
Gene Symbol	IL13ra1
Uniprot ID	I13R1_MOUSE
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081